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Q&A集

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2021/9/28

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The screenshot displays the Veritas software interface. The 'Reports' menu is highlighted in the top navigation bar. A dropdown menu is open, showing options: 'Antibody Custom', 'Antibody Screening/Identification', 'Antibody Screening Results', and 'Single Antigen Report'. The 'Antibody Custom' option is selected. Below this, the 'Custom Antibody Screening Report Setup' dialog box is open, titled 'HLA Fusion™'. The dialog box contains various sections for configuring the report, including 'Patient Information', 'Test Configuration', 'Tail Analysis Results', 'Overall Results/Assignments', 'Epitope Analysis Results', 'Test and Catalog Details', and 'Sample Information'. Each section has a 'Check All' and 'Uncheck All' button. The 'Setup' button in the bottom right corner of the dialog box is highlighted with a red box. An orange callout box at the bottom right of the dialog box contains the text: 'レポートに出力する項目を選択' (Select items to output in the report).

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Epitope Analysis

Class I	Specificity	TP	FN	Mean (Baseline) of Positives
	B45	1	0	26434.08
	B76	1	0	24681.3
	B44	2	0	23915
	A33	2	0	21937.27
	A68			
	A34			
	A26	1	0	19445.58
	A31	1	0	19219.48
	A66	2	0	18255.89
	A25	1	0	17949.33
	B82	1	0	17663.97
	A69	1	0	17441.25
	A74	1	0	16800.00
	A29			
	A3			

CREG Table

○ Red Circle: DSA Positive
○ Blue Circle: Non DSA Positive
○ Yellow Circle: DSA Negative

A Locus:
█ Strong Cross-reactivity
█ Cross-reactivity
 CREG (Cross-reactivity Group)

B Locus:
█ Strong Cross-reactivity
█ Cross-reactivity
 CREG (Cross-reactivity Group)

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Test Configuration

User Settings:
 Excluded Antigens, Positive Re-Threshold, Normalization Formula Level
 NC and PC Values, % PRA and Cutoffs by Region

Overall Results/Assignments

Overall Pos/Neg Other
 % PRA % Donor
 Antibody New
 Possible Epitope Amin

Epitope Analysis Results

Antibody Specificity (REQUIRED for this group)
 CREG Analysis Details:
 TP, FN, Mean (Raw) of Positives
 CREG Chart, Circle DSA

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×

Tail Analysis Results

Antibody Specificity (REQUIRED for this)

Tail Analysis Details:
TP/FP/FN/TN, R Value, Average Score, % Inclusion, Strength Index

Manual Tail Assignment

Test and Catalog Details

Bead ID or Test Well (REQUIRED for this)

Test Values: Rxn, Raw, Normalized, Count

SFI Raw, SFI Normal

Sero Specificity

Allele Specificity

Sort by: Bead ID ASC DESC

Test and Catalog Details

Bead ID or Test Well (REQUIRED for this)

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×

Tail Analysis Results

Antibody Specificity (REQUIRED for this)

Tail Analysis Details:

TP/FP/FN/TN, R Value, Average Score, % Inclusion, Strength Index

Manual Tail Assignment

Test and Catalog Details



Bead ID or Test Well (REQUIRED for this)

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