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FilmQA Pro - Case.2014-11-19_48bit - Calibration tool

File Panel Data Help

- Film calibration (ordinary)
 - Data - Calibration film
 - Tool - Calibration tool
- Dose map (single scan)
 - Data - dose film
 - Data - dose map (read only)
 - Data - consistency map (read only)
 - Dependence - Film calibration (ordinary)
- Dose to plan comparison
 - Data - Treatment plan
 - Data - Dose map (read only)
 - Data - Dose map registered (empty read or
 - Data - Comparison dose to plan (empty rea
 - Dependence - Dose map (single scan)
 - Tool - Image comparer

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The screenshot displays the FilmQA Pro software interface. The main window shows a dose map with a grid overlay. Two blue dashed boxes highlight specific areas on the map. The left sidebar contains a tree view of project data, with 'Dose map (single scan)' selected. The right sidebar shows a detailed information panel for the selected 'dose film'.

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The screenshot displays the FilmQA Pro software interface. The main window shows a grayscale image of a treatment plan with a grid overlay. A green dashed box highlights a specific region, and a white circle indicates the isocenter. The axes are labeled with numerical values: X-axis (0 to 280) and Y-axis (0 to 342).

The right-hand panel provides detailed information about the treatment plan:

- Info:** Name: Treatment plan; Dimensions: (X=left, Y=top) (W=width, H=height) □=Area; (X=0, Y=0) (W=280, H=94) □=26320 (pix); (X=0.0, Y=0.0) (W=547.7, H=282.0) □=1544.5 (mm.cm²); Size: 157.9 kB; Pixel format: 48 bpp RGB; Resolution: (X=13 dpi, Y=8 dpi); Light path: Undefined.
- Statistics:** Id: cc0e5e79-ff98-4b26-b5f2-2de9b7cd8619; File: L*\EBT3*Case#3; Origin: 愛知県がん-HP*愛知県がん 2014/11/19 16:41:30; Type: Plan.
- Profile:** Red quantity: L*\EBT3*Case#3; Green quantity: L*\EBT3*Case#3; Blue quantity: L*\EBT3*Case#3.
- Contour:** Channel range 'Red': [1.2, 160.0] cGy; Channel range 'Green': [1.2, 160.0] cGy; Channel range 'Blue': [1.2, 160.0] cGy; Source: Not available.
- Iso-map:** Frames: # S (X=left, Y=top) (W=width, H=height) □=Area; pix: 1 R (X=101, Y=29) (W=51, H=33) □=1704 mm.cm²; 1 R (X=198.0, Y=88.3) (W=100.0, H=100.0) □=100.0.
- Surface:** Fiducials: # (X=middle, Y=center) (W=width, H=height) □=Area Type; pix: 1 (X=140, Y=47) isocenter point mm.cm²; 1 (X=272.9, Y=139.5) isocenter point.
- Report:** Registration template: Not available.

The bottom status bar indicates: Image 'Treatment plan' - 185 %.

14:22> Panel switched to 'Image viewer and editor'.

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The screenshot shows the FilmQA Pro interface with a treatment plan loaded. The main window displays a grid with a central isocenter point. The left sidebar shows a tree view with 'Data - Treatment plan' selected. The right sidebar shows various tool icons and a detailed information panel. The information panel includes dimensions, fiducial data, and channel ranges.

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Patient Information:
DOB: Oct 21, 1941
Sex: M
ID: 0014033085
Oncologist: UNKNOWN, UNKNOWN
Disease: 12520

Plan Information:
Plan Label: Plan_01
Plan Status: Approved
Plan Date: Nov 6, 2014 6:23:14 PM
Position: HFS
DQA Plan: Plan_01

Delivery QA Plan Complete:
The procedures from this plan may be used to treat the phantom.

Delivery QA Setup / Analysis:
DQA Plan: Plan_01
Phantom Selector: 0210075
Retrieves Original Phantom
Change Phantom Couch
Laser Control:
 View Lasers
 Move Lasers
Save Position

Phantom Tools:
Restore Position
Move Phantom

Dose Calculation:
Calculation Grid: Fine
 Save Dose Per Angle
Start
Cancel

Displayed Image:
 Patient Phantom
Isodose Selector:
 Patient Plan Local
Restore From Plan
View Phantom Density Image
DQA Plan IVDT Selector

Dose Display:
 Isodose
64.2
60
57
54
48
42
30
18

ROI Display:
Tumor Settings

Name	Display	Color
PTV-60Gv	<input type="checkbox"/>	Magenta



Sensitive Structure Settings

Name	Display	Color
EAR INN R PRV	<input type="checkbox"/>	Cyan
EAR INN L PRV	<input type="checkbox"/>	Orange
LENS L	<input type="checkbox"/>	Grey
EYE L	<input type="checkbox"/>	Red
CTV-60Gv	<input type="checkbox"/>	Red
SPINAL_CORD	<input type="checkbox"/>	Green
BRAIN STEM	<input type="checkbox"/>	Cyan
SPINAL_CORD PRV	<input type="checkbox"/>	Green

Views:
- Transverse: 65, HFS
- Coronal: 93, HFS
- Sagittal: 121, HFS

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



The screenshot displays the FilmQA Pro software interface, specifically the 'Image comparer' window. The interface is divided into several sections:

- Left Panel (File, Panel, Data, Help):** A tree view showing project data including 'Dose map (single scan)', 'Data - dose film', 'Data - dose map (read only)', 'Data - consistency map (read only)', 'Dependence - Film calibration (ordinary)', 'Dose to plan comparison', 'Data - Treatment plan', 'Data - Dose map (read only)', 'Data - Dose map registered (empty read only)', 'Data - Comparison dose to plan (empty read only)', 'Dependence - Dose map (single scan)', 'Tool - Image comparer', 'Case object management', 'Add new case object', and 'New treatment case'.
- Central Image Area:** Shows two side-by-side images of a film scan. The left image is the 'Dose map (projected) Red' and the right is the 'Treatment plan Red'. A vertical axis on the left indicates distance in mm from 0.0 to 112.1. A horizontal axis at the bottom indicates distance in mm from 0.0 to 220.0. A path profile chart is visible on the right side of the image area.
- Registration Settings (Bottom Center):** A control panel for registration with the following values:
 - offset - X / mm: 82.0
 - offset - Y / mm: 125.2
 - rotation angle: 0.4 °
 - max reference: 160.0 cGy
 - Other offset and rotation fields are set to 0.0 mm.
- Image Comparison Settings (Right Side):** A control panel for image comparison with the following values:
 - Distance: 3.0 mm
 - Threshold: 100.0 %
 - Default value: 100.0 %
 - Passing rate 'red': 82.86 %
 - Tolerance: 3.0 mm
 - Distance: 3.0 mm
 - Threshold: 100.0 %
- Bottom Status Bar:** Displays the registration details: '14:07> Registrant projected at {W=280, H=94} using coordinate system 'd'={X=74.72015, Y=-4.544569}, a=0.2399948°, w={17.31 %, 11.29 %}' and 'Bi-linear' interpolation.'

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








Single Channel Correction (red), 3 mm/3%, TH10%, Area 15x15 cm²

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