

## 【正誤表】

P12の”表 1. EBT3 の組成”を以下の表と差し替える

## 【誤】

Material	Thickness	Approx. imate density	COMPOSITION (ATOM%)										
	microns	g/cm <sup>3</sup>	H	Li	C	N	O	S	Ba	Bi	Na	Cl	Al
Clear Polyester film base	125	1.35	36.40%	0.00%	45.50%	0.00%	18.20%						
Active layer (contains approx. imately 5% moisture)	28	1.35	55.80%	0.60%	31.10%	0.40%	10.50%	0.10%			0.10%	0.10%	1.40%

## 【正】

Material	Thickness	Approx. density	COMPOSITION (ATOM%)										
	microns	g/cm <sup>3</sup>	H [%]	Li [%]	C [%]	N [%]	O [%]	S [%]	Ba [%]	Bi [%]	Na [%]	Cl [%]	Al [%]
Clear polyester film base	125	1.35	36.4		45.5		18.2						
Active layer (contains approx. imately 5% moisture)	28	1.20	56.9	0.6	27.6		13.3						1.6
Total	278	1.35	38.4	0.1	43.7		17.7						0.2