


# One Lambda HPへのアクセス

- <https://www.thermofisher.com/onelambda/wo/en/home.html>にアクセス
- Loginをクリック



- 登録したEメールアドレスとパスワードを入力し、Sign Inをクリック




The screenshot shows a 'Sign In' form with the following elements:

- Sign In** (Title)
- Instruction: Please sign in to access your account. Customers must log in to access some content, recorded webinars, training materials, documents, and software.
- Work Email \*** label above a text input field. An orange box with the text **Eメールアドレス** is overlaid on the input field.
- Password \*** label above a text input field. An orange box with the text **パスワード** is overlaid on the input field.
- [Forgot Password / Reset Password](#) link.
- Cancel** button (dark grey).
- Sign In** button (maroon). A yellow arrow points to this button.
- Footer text: Need Help? Call website support Toll Free: +1 800 822 8824 or [Email support team](#).
- Bottom banner: Don't have a profile? **Register Now**.


# ダウンロードページへのアクセス

- Customer Care をクリック
- Documentation&Software をクリック



The screenshot shows the ONE LAMBDA website header. The logo on the left reads "ONE LAMBDA A Thermo Fisher Scientific Brand". On the right, there is a user greeting "Welcome, [redacted] Log Out" and a navigation menu with "Products & Services", "Customer Care", "Learning Center", and "Company". A search icon is also present. The "Customer Care" menu is open, displaying a list of options: "Products & Support Files", "Technical Product Support", "Documentation & Software" (highlighted with a red box), "Knowledge Base Articles", and "RSS Feeds". Each item has a lock icon to its right. Below the navigation, a banner for "Essary of the One 3Screen Single Antigen Bead" is visible.

- Product Lineからダウンロードしたい資料の対象製品を選択




The image illustrates the process of selecting a target product from a list. On the left, a dropdown menu is open for 'Product Line', showing a list of products. The 'EN-English' option is selected. An orange arrow points to the 'Product Line' dropdown. On the right, the expanded dropdown menu for 'LABType SSO' is shown, with a yellow hand icon pointing to the 'LABType SSO' option. An orange arrow points from the left panel to the right panel.

| Product Line   | EN-English | Type | Lot/Version |
|--|------------|------|-------------|
| <input type="checkbox"/> Adsorb Out                      |            |      |             |
| <input type="checkbox"/> AllSet Gold SSP High-Resolution |            |      |             |
| <input type="checkbox"/> AllSet Gold SSP Low-Resolution  |            |      |             |
| <input type="checkbox"/> AllType FASTplex                |            |      | Ce          |
| <input type="checkbox"/> AllType NGS Ancillary Reagents  |            |      |             |
| <input type="checkbox"/> AllType NGS Assays              |            |      | Ce          |
| <input type="checkbox"/> AmbiStrips                      |            |      |             |
| <input type="checkbox"/> Ancillary Products              |            |      |             |
| <input type="checkbox"/> Ancillary Reagents              |            |      | Ce          |


| LABType SSO   | EN-English | Type | Lot/Version |
|---|------------|------|-------------|
| <input type="checkbox"/> LABType Ancillary Reagents |            |      |             |
| <input checked="" type="checkbox"/> LABType SSO     |            |      |             |
| <input type="checkbox"/> LABType SSO HD             |            |      |             |
| <input type="checkbox"/> LABType XR and CWD         |            |      |             |
| <input type="checkbox"/> LABXpress Pipettor         |            |      |             |
| <input type="checkbox"/> Lambda Monoclonal Trays    |            |      |             |
| <input type="checkbox"/> Lambda Scan Software       |            |      |             |
| <input type="checkbox"/> LAT                        |            |      |             |
| <input type="checkbox"/> LAT Mixed                  |            |      |             |

# 対象ファイルの選択


- Typeからダウンロードしたい資料を選択




テンプレートファイルの場合：  
Luminex Templatesを選択



カタログファイルの場合：  
Software Analysis Filesを選択




ワークシートの場合：  
Worksheets/Datasheetsを選択



# 対象ロットの選択

- Lot/Versionから製品のロットを選択



- Downloadボタンをクリックしてダウンロード

| Title   | Type                    | Publish Date | Catalog ID     | Lot/Version | Download  |
|---|-------------------------|--------------|----------------|-------------|---|
| LABType SSO Class I A Locus Typing Test, Lot 018 - HLA Catalog File | Software Analysis Files | 05/15/2023   | RSSO1A, RSO1AT | 018         |  |

